

BED SORE PREVENTION PAD (BIG SQUARE)



ABOUT THE PRODUCT

A bed sore prevention pad with big square tiles is designed to reduce pressure and protect the skin from ulcers. The larger gel or foam squares offer broader support, helping distribute body weight evenly and reducing stress on high-risk areas like the hips, tailbone, and heels. Its open-grid design promotes airflow, keeping the skin cool and dry to prevent irritation. This pad works well on beds, wheelchairs, and chairs for effective, comfortable pressure relief.

BENEFITS

- Small square tile gel layer 10 mm thick with cloth cover provides ultimate comfort and better posture.
- Reduces Skin temperature, heat flux and relative humidity.
- Temperature remains static under the ischial tuberosities and posterior thighs.
- Effectively reduce ischemic tissue damage.
- Reduces the incidence of pressure sores.

HOW TO USE

- Place the pad under the pressure-prone area (heels, elbows, hips, tailbone).
- Make sure it lies flat with the cushioned side facing up.
- Position the patient gently so the area rests comfortably on the pad.

Instructions : If the product gets dirty or sticky, clean with talcum powder, rinse with mild shampoo and water, shade dry, and avoid direct heat.

Safety : Keep away from heat and sharp objects.

THANK YOU FOR YOUR PURCHASE — YOUR COMFORT MATTERS TO US.

**FEEL FREE TO CONTACT US FOR MORE INFORMATION
WWW.LIVISH.IN
PH NO- +91 99620 94754**