

WHEEL CHAIR SUPPORT (ECONOMY)



ABOUT THE PRODUCT

A Wheel Chair Support labeled "1 Size Fits All" is a versatile cushion designed for universal application across most standard wheelchairs and seating. Typically constructed with a durable foam core or a simple gel/foam layer, its main advantage is ease of use and broad compatibility, aiming to provide a basic level of pressure relief and enhanced comfort. While it offers fundamental padding and improves weight distribution compared to a bare seat, its non-specialized sizing means it prioritizes simplicity and broad utility over highly customized anatomical fit.

BENEFITS

- 15 mm gel layer in a profiled PU foam base with cloth cover helps reduce pressure points for superior comfort and posture.
- It lowers skin temperature, heat, and humidity, keeping temperature stable under the hips and thighs to prevent tissue damage.
- Helps reduce the risk of pressure sores and stays durable, long-lasting, and infection-free.

HOW TO USE

- Position the cushion on the wheelchair seat with the cloth side facing up.
- Sit in the center of the cushion, ensuring your hips and thighs are evenly supported.
- Shift the cushion to align with your comfort and posture. Ensure it doesn't slip before use.

Instructions : If the product gets dirty or sticky, clean with talcum powder, rinse with mild shampoo and water, shade dry, and avoid direct heat.

Safety : Keep away from heat and sharp objects.

THANK YOU FOR YOUR PURCHASE — YOUR COMFORT MATTERS TO US.

**FEEL FREE TO CONTACT US FOR MORE INFORMATION
WWW.LIVISH.IN
PH NO- +91 99620 94754**