

WHEEL CHAIR SUPPORT (ROUND BALLS)



ABOUT THE PRODUCT

A wheelchair support pad featuring round gel balls is a sophisticated cushion designed for superior pressure relief and comfort for wheelchair users. This pad consists of numerous individual, round gel or fluid-filled balls encased in a flexible cover. This construction allows the pad to constantly conform and adjust to the user's changing contours and movements, distributing weight exceptionally well and eliminating pressure hot spots. The spaces between the balls promote continuous airflow, which helps manage heat and moisture, making it highly effective for preventing pressure sores during prolonged sitting.

BENEFITS

- 20 mm round gel layer with 1-inch soft PU foam and cloth cover for superior comfort & posture.
- Helps reduce skin temperature, heat buildup, and humidity.
- Keeps temperature stable under the hips and thighs.

HOW TO USE

- Place the cushion on the wheelchair seat with the round gel side facing up.
- Sit gently, ensuring your hips and thighs are centered on the cushion.
- Adjust your position for balanced pressure relief and comfort.

Instructions : If the product gets dirty or sticky, clean with talcum powder, rinse with mild shampoo and water, shade dry, and avoid direct heat.

Safety : Keep away from heat and sharp objects.

THANK YOU FOR YOUR PURCHASE — YOUR COMFORT MATTERS TO US.

**FEEL FREE TO CONTACT US FOR MORE INFORMATION
WWW.LIVISH.IN**

PH NO- +91 99620 94754