


# BED SORE PREVENTION PAD (SMALL SQUARE)



## ABOUT THE PRODUCT

A bed sore prevention pad with small square tiles is a therapeutic overlay designed to prevent pressure ulcers. Its grid of individual gel or foam squares allows the pad to mold to the body, evenly distributing weight and reducing pressure on bony areas like the tailbone, hips, and heels. The spaces between the tiles improve airflow, keeping skin cool and dry to prevent breakdown. This pad is ideal for wheelchairs, chairs, or targeted support on beds.

## BENEFITS

- Small square tile gel layer 10 mm thick with cloth cover provides ultimate comfort and better posture.
- Reduces Skin temperature, heat flux and relative humidity.
- Temperature remains static under the ischial tuberosities and posterior thighs.
- Effectively reduce ischemic tissue damage.
- Reduces the incidence of pressure sores.
- Last long, washable and infection free.

## HOW TO USE

- Place the pad under the pressure-prone area (heels, elbows, hips, tailbone).
- Make sure it lies flat with the cushioned side facing up.
- Position the patient gently so the area rests comfortably on the pad.

**Instructions :** If the product gets dirty or sticky, clean with talcum powder, rinse with mild shampoo and water, shade dry, and avoid direct heat.

**Safety :** Keep away from heat and sharp objects.

**THANK YOU FOR YOUR PURCHASE — YOUR COMFORT MATTERS TO US.**

**FEEL FREE TO CONTACT US FOR MORE INFORMATION  
WWW.LIVISH.IN  
PH NO- +91 99620 94754**