

# HEEL CUP



## ABOUT THE PRODUCT

Heel cup gels are a comforting solution for foot support that you can use at any time of the day. They are small, cushioned, and perfectly designed, making them ideal for relieving discomfort, reducing pressure, or just providing extra comfort with your shoes. With various styles and materials available, from classic silicone to advanced gel technologies, there is always a heel cup gel to suit your needs. Whether you choose them for daily wear or specific activities, a heel cup gel always brings relief to your feet.

## BENEFITS

- Help in correcting the alignment while walking, absorbs shock and give right elevation to the heel.
- Prevent feet from drying, cracking and peeling. Provides extreme cushioning to heel for relieving Heel Pain, Plantar Fasciitis & Achilles tendinitis.
- Medical grade elastomeric cup contoured to cradle heel cracks & ankle together, thereby prevents side-to-side motion stabilizing foot as it strikes on ground.

## HOW TO USE

- Place the heel cup at the back of your shoe
- Ensure it sits flat and secure
- Wear shoes and adjust for comfort
- Use in both shoes for balanced support

**Instructions :** If the product gets dirty or sticky, clean with talcum powder, rinse with mild shampoo and water, shade dry, and avoid direct heat.

**Safety :** Keep away from heat and sharp objects.

**THANK YOU FOR YOUR PURCHASE — YOUR COMFORT MATTERS TO US.**

**FEEL FREE TO CONTACT US FOR MORE INFORMATION  
WWW.LIVISH.IN**

**PH NO- +91 99620 94754**